

3-2016

Surface functionalized nylon capillary-channeled polymer (C-CP) fibers for protein ion-exchange separations

Liuwei Jiang
Clemson University

R. Kenneth Marcus
Clemson University

Follow this and additional works at: <https://tigerprints.clemson.edu/cars>

 Part of the [Chemistry Commons](#)

Recommended Citation

Jiang, Liuwei and Marcus, R. Kenneth, "Surface functionalized nylon capillary-channeled polymer (C-CP) fibers for protein ion-exchange separations" (2016). *Chemistry Annual Research Symposium*. 14.
<https://tigerprints.clemson.edu/cars/14>

This Poster is brought to you for free and open access by the Student Works at TigerPrints. It has been accepted for inclusion in Chemistry Annual Research Symposium by an authorized administrator of TigerPrints. For more information, please contact kokeefe@clemson.edu.



Surface functionalized nylon capillary-channeled polymer (C-CP) fibers for protein ion-exchange separations

Liuwei Jiang and R. Kenneth Marcus

CLEMSON UNIVERSITY, Department of Chemistry, Clemson, SC 29634



What are C-CP fiber columns ?




Figure 1. C-CP fiber column construction process

C-CP fibers are made from melt extrusion of commonly used polymers (e. g. PP, PET and nylon). There are 8 capillary channels on C-CP fiber which greatly increases the surface area and fluid transportation. C-CP columns are packed by pulling fibers through PEEK tubing.

Virtues of C-CP fiber column:


- Highly permeable (a 200 mm×2.1 mm i.d. C-CP column can be run at **10 mL min⁻¹** flow rate with back-pressure <1500 psi)
- Fast protein purification at high solvent flow rates
- **Low cost**




C-CP fibers on spool

C-CP fiber column

Cross-section of column




C-CP fibers under SEM



Cross-section of C-CP fibers under SEM

Nylon C-CP fiber Surface Modification

Native nylon C-CP fibers have a very low density (~20 μmol g⁻¹) of carboxylic acid groups that limits its performance on protein separations. To overcome this challenge, surface modification was done on native nylon fibers by grafting monomer ligands (acrylics) via microwave-assisted free radical polymerization.



- **Simple modification:** Modification is fast (~10 min) and low-cost using a domestic microwave oven
- **High ligand density of modified fibers:** -COOH or -SO₃H ligands were grafted on native C-CP fibers in up to 500 μmol g⁻¹ density
- **Versatile modification:** Different types of acrylic monomers are industrially available.

Protein Separations on C-CP fibers




Figure 2. Protein separations comparison on nylon-C-CP fiber column before and after surface modification

- Proteins: 1. myoglobin, 2. α-chymotrypsinogen A, 3. cytochrome C, 4. lysozyme
- Column size: 200 mm length × 0.762 mm.i.d.
- Mobile phase A: 20 mM phosphate pH 6.5, Mobile phase B: 1.0 M NaCl in A
- Flow rate: 0.1 mL min⁻¹ (Black line), 0.5 mL min⁻¹ (Red line)

Before surface modification

- Very poor protein ion-exchange separation on native fibers
- Co-elution of proteins

After surface modification

- **(Black line):** Improved ion-exchange separation on modified column (Flow rate: 0.1 mL min⁻¹)
- **(Red line):** Ion-exchange separation at a 5× increase flow rate (Flow rate: 0.5 mL min⁻¹)
- Increase of flow rate by 5× does NOT impair the separation resolution
- Decrease of total separation time from 30 min to 6 min

Protein Loading/Elution

Table 1. Dynamic binding capacity (DBC) of lysozyme on nylon C-CP fiber columns at different protein loading concentrations.

Loading Concentration (mg mL ⁻¹)	Native Nylon DBC (mg g ⁻¹)	Modified Nylon DBC (mg g ⁻¹)
0.05	0.63±0.44	19.13±0.18
0.1	1.81±0.02	20.11±0.02
0.2	0.87±0.19	20.56±0.01
0.4	0.81±0.27	21.07±0.38
0.6	1.24±0.47	22.22±0.03
0.8	0.94±0.67	22.28±0.04
1.0	1.66±0.10	22.60±0.89

Figure 3. Consistency of lysozyme loading and elution

- 10 lysozyme loading and elution chromatograms. No column regeneration was carried out between runs
- Excellent lysozyme recovery on modified nylon column
- **RSD of binding capacity = 0.3% (n=10)**





Figure 4. Column to column consistency

- Negligible variations in breakthrough curves
- Batch-to-batch reproducibility of the modification could be improved by using more controllable microwave instruments
- **RSD of binding capacity = 3% (n=5)**



Conclusions

Microwave-assisted grafting polymerization largely improved the performance of the nylon C-CP fiber packed column in protein separations. Applying this simple but versatile modification method on nylon C-CP fibers offers cost-efficient LC stationary phase that are capable of fast protein separations and high-throughput protein processing.

Acknowledgments

This material is based upon work supported by the National Science Foundation Division of Chemistry under Grant No. CHE-1307078.

References

Please contact liuwei@clmson.edu for a reference list.